

Wet Chemistry Consumables

Kjeldahl Catalysts and Glass Tubes

Increase your Kjeldahl Digestion Efficiency with Catalyst Tablets. Consult OPSIS LiquidLINE laboratory to select the right combination for your application method. Select among our portfolio of Borosilicate Glassware. All consumables are made in Europe.



KT-210-A
Missouri 5g
Cu 0.3%



KT-211-A
Missouri 3.5g
Cu 0.3%

Recommended



KT-220-A
Micro Copper
1.65g
Cu 9%



KT-230-A
Universal
Copper 5g
Cu 9%



KT-240-A
Titanium 3.71g
Cu-TiO₂



KT-250-A
Selenium 3.5g
Se 0.1%



KT-260-A
Wieninger 3.5g
Se-Cu 2%



TU-220-A
Conical Flask
250 ml
10 flasks



TU-214-A
Tubes 250 ml
20 tubes



TU-223-A
Tubes 450 ml
10 tubes



TU-221-A
Tubes 750 ml
4 tubes



TU-222-A
Boiling Rods
10 rods

Extraction Thimbles and OPSIS LiquidLINE Filters

Select among our portfolio of standard extraction thimbles, compatible with all instruments on the market. Use our unique filters to increase efficiency of your Hydrolysis and Fiber methods.



SX-040-A
Cellulose
Thimble 33x80
25 pcs



SX-041-A
Cellulose
Thimble 25x80
25 pcs



SX-044-A
Cellulose
Thimble 40x80
25 pcs



SX-045-A
Glass Microfiber
Thimble 33x80
25 pcs



SX-046-A
Glass Microfiber
Thimble 40x80
25 pcs



SX-043-A
HydROC Filters
50 filters



SX-050-A
Crude Fiber
Bags 25/17
25µm porosity
50 pcs

OPSIS LiquidLINE Accessories

Accessories to your KjelROC Distillation and Analyzer System

Configure your KjelROC system with additional accessories to lower risk for mistakes, increase precision and for added efficiency.



LA-020-A
LabConnect
LIMS software
to increase
efficiency



KD-002-A
Recovery Test
with delivery of
unit



KD-026-A
IQ/OQ/PQ
Test Documents



KD-020-A
Flexible Tank
Level Sensors
4 sensors



KD-022-A
Burette Kit for
quick exchange
of Titrant
solution



KD-023-A
Büchi Adapter
for using Büchi
glass tubes in
Analyzer



KD-025-A
Tube Adapter
for using 100 ml
glass tubes in
Analyzer



KD-024-A1/2/3
External Titrator
Package
1 cable



KD-030-A1
Glass Splash-
head for
Standard
Analysis



KD-030-A3
PP Splashhead
for Analyzer and
High Volume
analysis



KD-001-A1
UK Power
Cable

Accessories to your KjelROC Digester


Increase your Kjeldahl Digestion efficiency with added accessories.




DI-001/002-A
Tube Rack 10
or 20 position
Digester



DI-030-A1/A2
Exhaust 10
or 20 position
Digester



DI-050-A1/A2
Fixed cooling
Stand for 10
or 20 position
Digester



KD-051-A
Cooling Stand



DI-052-A2
Washing Lid
for 20 position
Digester



DI-054-A1/A2
Water Jet Pump
for 10 or 20
position
Digester



DI-111-A
Water Trap for
Liquid Digestion
and Scrubber

LP 1024
2021 07