

# OPSIS LiquidLINE KjelROC Scrubber

Environment-friendly Digestion



# KJELROC SCRUBBER

A Scrubber is recommended for fume control and in areas where water is expensive. The KjelROC Scrubber can save costs.

## RELIABLE AND EASY WITH A SMALL FOOTPRINT

The KjelROC Scrubbers are built to work in an environment with an excess of corrosive fumes and high heat. All materials are carefully selected.

- Cooler and containers in Borosilicate glass and connection tubes in Viton, in order to withstand the demanding environment and reduce maintenance costs
- Small footprint makes it easy to install the Scrubber close to the KjelROC Digestor
- Easy access to exchange containers, which are all made in high quality Borosilicate Glass.
- Accessible 6-step suction control to regulate suction power



## NO WATER IS NEEDED

The KjelROC Scrubber does not require any tap water and there is a build-in pump for suction power.

## OPTIONS

All KjelROC Scrubbers can be supplied with optional configuration

- Water cooler, Water Accessory to even further increase the efficiency of the fume extraction. Can also be used with activated carbon filter.
- Water Trap for digestion of liquid samples

Opsis AB, Box 244  
SE-244 02 Furulund Sweden  
Telephone +46 46 72 25 00  
Telefax +46 46 72 25 01  
E-mail [info@opsis.se](mailto:info@opsis.se)  
[www.liquidline.se](http://www.liquidline.se)



LP 1018  
2021 02